



# ULTRA DISINFECTANT

Virucidal activity against enveloped viruses including Coronavirus (BCoV) (as a surrogate for other members of the Coronavirus family including MERS-CoV).



High performance disinfectant with broad spectrum activity against bacteria, yeast and viruses.

## KEY FEATURES:

- Aldehyde-free disinfectant for use on most washable hard surfaces.
- Independently tested to European Standard Test Method EN 1276:2019, EN 14476 and EN 13697.
- Based on dodecyldimethylammonium chloride (DDAC).
- Suitable for the disinfection of hard surfaces in hospitals, long term care facilities, institutions, restaurants, food processing plants and schools.

High performance disinfectant



Refill trigger spray bottles from 5 litre bottle to reduce plastic waste and achieve better economy in use.



## HIGH EFFICACY:

- Bactericidal (inc. MRSA, E. coli)
- Yeasticidal
- Virucidal (against enveloped & non-enveloped viruses including, Influenza and Bovine Corona Virus (BCoV) (as a surrogate for other members of the Coronavirus family including MERS-CoV).

Claim	European Standard Test Method	Contact time
Bactericidal	EN 1276:2019 (inc. MRSA, E. coli)	30 seconds
	EN 13697 (inc. E. coli)	5 min
Yeasticidal	EN 16615	1 min
Virucidal	EN 14476 (inc. Influenza)	5 min
	EN 14476 (inc. Coronavirus)	1 min