

# MEDIANA HEARTON AED A15



## SPECIFICATIONS

### DEFIBRILLATION ELECTRIC SHOCK

Waveform:	Biphasic Truncated Exponential (BTE) waveform (Impedance compensation)
Energy:	Adult: 185 to 200J ( $\pm 5\%$ ) Paediatric: 45 to 50J ( $\pm 5\%$ )
Operating Mode:	Semi-Auto
Charge Time:	Less than 12 seconds
Analysis Time:	8-13 seconds
Shock Time:	Delivers 200 shocks with full battery

### ECG

Lead:	II (RA, LL)
Patient Impedance:	25 to 175ohm
Detection V/F	more than 200 $\mu$ V
V/T	more than 160 per min

### INDICATION

#### CONTROLS

Standard:	Open (power) switch, shock button, Adult/Paediatric mode switch
-----------	---

#### INDICATORS

Visible:	ICON Indicator, Status LCD (unit status, battery status, temperature status), LED (Adult/Paediatric mode switch LED)
Audible:	Speaker (voice prompt), Beep (CPR indication)

### PHYSICAL

Dimensions:	294H x 240W x 95mmD
Weight:	Approx. 2.65kg (including battery, excluding pads)

### SELF TEST

Cycle:	Every 24 hours, 1 week, 1 month Power on self test, battery insertion self test
Test Result:	Status LCD displays 'O' / 'X'

### DATA BACKUP AND COMMUNICATION

Standard SD card, IR communication port
---

### ACCESSORY SPECIFICATIONS

#### ADULT/PAEDIATRIC PADS

Standby Life:	2 years from date of manufacture
Electrodes:	Disposable pads
Placement:	Adult: Anterior-lateral Paediatric: Anterior-posterior
Minimum active gel area:	80cm <sup>2</sup> +/-5%
Cable Length:	Approx. 1.8m

#### Environmental Conditions

Temperature:	Operation: 0 to 43°C (32 to 109.4°F) Storage: 0 to 43°C (32 to 109.4°F)
Relative Humidity:	5 to 95% RH (non-condensing)

#### BATTERY

Type:	LiMnO <sub>2</sub> , disposable, long-life primary cell
Voltage/Capacity:	15V, 4200 mAh
Shelf Life:	2 years from date of manufacture (in the original packaging)
Standby Life:	5 years from date of manufacture (inserted in the AED)
Discharge:	A minimum of 200 shocks (except the CPR period between the defibrillation therapy) or 10 hours of operating time at 20°C

#### Environmental Conditions

Temperature:	Operation: 0 to 43°C (32 to 109.4°F) Storage: 0 to 43°C (32 to 109.4°F)
Relative Humidity:	5 to 95% RH (non-condensing)

Information correct at time of going to press. Subject to change without notice.

[www.medianadefib.co.uk](http://www.medianadefib.co.uk)

t. +44 (0) 08456 448808  
e. [sales@reliancemedical.co.uk](mailto:sales@reliancemedical.co.uk)  
w. [www.medianadefib.co.uk](http://www.medianadefib.co.uk)

